


# PureFlex<sup>®</sup> Arterial Cannulae



A complete range of adult arterial cannulae crafted to meet diverse clinical needs in conventional cardiac surgery

## New, Curved and Straight Tip Arterial Cannulae

The all new family of curved and straight tip PureFlex Arterial Cannulae have been designed to meet diverse clinical needs with different configurations. These cannulae provide smooth no-jet flow and a low pressure drop\*, while the new tip design grants ease-of-insertion. The upgraded, sleek DEHP-free tubing is highly transparent for greater visibility and its new wire-reinforced, flexible spring design, with reduced pitch, allows no kinking.



New, highly transparent  
DEHP-free tubing

New, flexible  
spring design  
allows no kinking

New, highly  
transparent tip  
design grants  
ease-of-insertion



## Fits all Adult Patient Sizes

The PureFlex Arterial Cannulae comes in 4 sizes and 36 different configurations with a 3/8" tubing connection, designed to fit all adult patients. The transparent, straight and curved tips come with either a suture flange or collar to accommodate all types of insertion procedures, while the more generous clamping area assures easier clamping during priming.



**PureFlex Arterial Cannulae**  
18, 20, 22, 24 Fr sizes available


## Pressure Drop\*



\* All data extracted from relative Instructions For Use

## Three Different Connections; Greater Versatility

Without luer lock · With luer lock for venting · With porous vent plug and no connector



### Available Configurations

#### Bent tip with collar

RAXX-1BB



Bent tip with collar, with connector 3/8" with luer lock

RAXX-1BC



Bent tip with collar, with connector 3/8" without luer lock

RAXX-1BP



Bent tip with collar, with 3/8" connection and vent plug

#### Bent tip with suture flange

RAXX-2BB




Bent tip with suture flange, with connector 3/8" with luer lock

RAXX-2BC



Bent tip with suture flange, with connector 3/8" without luer lock

RAXX-2BP



Bent tip with suture flange, with 3/8" connection and vent plug

#### Straight tip with collar

RAXX-1SB



Straight tip with collar, with connector 3/8" with luer lock

RAXX-1SC



Straight tip with collar, with connector 3/8" without luer lock

RAXX-1SP



Straight tip with collar, with 3/8" connection and vent plug

# LivaNova

Health innovation that matters

Discover more at [cannulae.livanova.com](http://cannulae.livanova.com)

[www.livanova.com](http://www.livanova.com)



FDA Cleared

Please always refer to the Instructions For Use (IFU) manual provided with each product for detailed information, warnings, precautions and possible adverse side effects.

Manufactured by:

**Sorin Group Italia Srl**

A wholly-owned subsidiary of LivaNova PLC  
Via Statale 12 Nord, 86 - 41037 Mirandola (MO) Italy  
Tel: +39 0535 29811 - Fax: +39 0535 25229  
[info.cardiacsurgery@livanova.com](mailto:info.cardiacsurgery@livanova.com)

For professional use. Please contact us through our website to receive instructions for use containing full prescribing information, including indications, contraindications, warnings, precautions and adverse events. Not approved in all geographies. Consult your labeling.

© 2019 LivaNova all rights reserved.