


## Arterial Cannulae, 135° Curved Tip, with Suture Collar, Wire-reinforced Tubing, 12 Fr/3.8 mm, A272-40N



Livanova 135° Curved Tip Arterial Cannulae allow for a high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	A272-40N
SIZE	12 Fr/3.8 mm
CONNECTION SITE	3/16"-1/4"
EFFECTIVE LENGTH (CM)	22
FEATURE	Without Connector
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10