

Suction Products, Suction Wands, Rigid Fluted Tip with Plastic Handle, 21 Fr, S900-17



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|------------------|-----------------------------------|
| ORDER CODE | S900-17 |
| SIZE | 21 Fr |
| BARBED CONNECTOR | 1/4" |
| AVAILABLE COATED | Depending on country registration |
| QTY PER BOX | 10 |
| NOTE | 2 pieces of S351-70 per package |