


## Thoracic Catheters, Straight, 12 Fr, TC-10112



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-10112
SIZE	12 Fr
CONNECTOR ACCEPTANCE	3/16 in
EFFECTIVE LENGTH (CM)	51,8
DRAINAGE HOLES	4
LENGTH OF PERFORATION	2.0 in / 5.0 cm