

## Thoracic Catheters, Right Angle, 16 Fr, TC-11116



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-11116
SIZE	16 Fr
CONNECTOR ACCEPTANCE	1/4 in
EFFECTIVE LENGTH (CM)	51,8
DRAINAGE HOLES	4
LENGTH OF PERFORATION	2.5 in / 6.0 cm