

## Thoracic Catheters, Right Angle, 24 Fr, TC-11124



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-11124
SIZE	24 Fr
CONNECTOR ACCEPTANCE	1/4 in
EFFECTIVE LENGTH (CM)	51,8
DRAINAGE HOLES	5
LENGTH OF PERFORATION	3.0 in / 7.5 cm