

## Thoracic Catheters, Right Angle, 28 Fr, TC-11128



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-11128
SIZE	28 Fr
CONNECTOR ACCEPTANCE	3/8 in
EFFECTIVE LENGTH (CM)	51,8
DRAINAGE HOLES	6
LENGTH OF PERFORATION	4.5 in / 11.0 cm