

Thoracic Catheters, Right Angle, 32 Fr, TC-11132



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2" increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-11132
SIZE	32 Fr
CONNECTOR ACCEPTANCE	3/8 in
EFFECTIVE LENGTH (CM)	51,8
DRAINAGE HOLES	6
LENGTH OF PERFORATION	4.7 in / 12.0 cm