


Venous Return Cannulae, Single Stage, Bullet Tip Venous Return Cannulae, Wire-reinforced Tubing, 16 Fr, V122-16



Livanova Pediatric Venous Return Cannulae are intended for venous direct cannulation of pediatric patients and are available with a single stage tip design. Tubings are designed in clear and transparent PVC which maintains elasticity during hypothermia and optimize flow characteristics. They feature a 90° angled tip which ensures optimal venous drainage. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|-----------------------|-----------------------------------|
| ORDER CODE | V122-16 |
| SIZE | 16 Fr |
| CONNECTOR ACCEPTANCE | 1/4" |
| EFFECTIVE LENGTH (CM) | 30 |
| AVAILABLE COATED | Depending on country registration |
| QTY PER BOX | 10 |