


## Venous Return Cannulae, Single Stage, Bullet Tip Venous Return Cannulae, Wire-reinforced Tubing, 20 Fr, V122-20



Livanova Pediatric Venous Return Cannulae are intended for venous direct cannulation of pediatric patients and are available with a single stage tip design. Tubings are designed in clear and transparent PVC which maintains elasticity during hypothermia and optimize flow characteristics. They feature a 90° angled tip which ensures optimal venous drainage. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	V122-20
SIZE	20 Fr
CONNECTOR ACCEPTANCE	1/4"
EFFECTIVE LENGTH (CM)	30
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10