


MICS Cannulae, Easy Flow DUO Femoral Arterial, 23 Fr, 103-310



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

ORDER CODE	103-310
SIZE	23 Fr
OUTER DIAMETER (MM)	7.7
TIP	Dispersion Flow
CONNECTION SITE	3/8"
EFFECTIVE LENGTH (CM)	28,5
AVAILABLE COATED	no
QTY PER BOX	5