


## MICS Cannulae, Vascular Dilator Kit, 8,12,16,20, 24 Fr, 200-120



Livanova Vascular Dilator Kit, for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035 x 180 cm guide-wire, a blade scalpel and 8,12,16,20 and 24 Fr dilators.

ORDER CODE	200-120
SIZE	8,12,16,20, 24 Fr
SCALPEL	yes
GUIDEWIRE	0.035 (in) x 180 cm
NEEDLE	18 Ga
QTY PER BOX	1