


## MICS Cannulae, RAP Femoral Venous, 23 Fr Distal/25 Fr Proximal, 200-150



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

ORDER CODE	200-150
SIZE	23 Fr Distal/25 Fr Proximal
OUTER DIAMETER (MM)	7.7/8.3
CONNECTION SITE	3/8" and 1/2"
EFFECTIVE LENGTH (CM)	65,5
AVAILABLE COATED	no
QTY PER BOX	5