


## MICS Cannulae, FlexFlow Venous, 23 Fr, 200-200



Livanova FlexFlow Venous Cannulae is a minimally invasive venous cannula with integrated obturator and malleable stylet designed for direct right atrium cannulation. The 6 cm holed area and the thin wire-reinforced tubing wall design provide high flow rates when using vacuum assisted venous drainage.

ORDER CODE	200-200
SIZE	23 Fr
OUTER DIAMETER (MM)	7.7
SIDE HOLES	yes
CONNECTION SITE	3/8"
EFFECTIVE LENGTH (CM)	28,5
AVAILABLE COATED	no
QTY PER BOX	5