


MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019



Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way.

Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

ORDER CODE	203-019
SIZE	19 Fr
OUTER DIAMETER (MM)	6.33
CONNECTION SITE	3/8"
EFFECTIVE LENGTH (CM)	26,7
AVAILABLE COATED	no
QTY PER BOX	5
NOTE	Kit Includes: Bi-Flow ECMO 19 Fr, Dilators (12, 19 Fr), Needle 18 Ga, Scalpel, Syringe 5 ml, Guidewire (100 cm, 0.035"), Cannula

