


## Advanced Circulatory Support, Extracorporeal Life Support , Bi-Flow ECMO, Femoral Cannula 19Fr, 303-019



BI-FLOW ECMO is the only bidirectional arterial cannula in the market designed to prevent leg ischemia in an easy, safe and reproducible way.

Bi-Flow ECMO obviates the need for separate downstream perfusion or other techniques to protect against lower limb ischemia

Bi-Flow ECMO is unique in the market . Its patented design incorporates a unique elbow and downstream perfusion channel.

Bi-Flow ECMO is PH.I.S.I.O. coated. The PH.I.S.I.O. coating improves blood compatibility of the device by reducing platelet adhesion to the coated surfaces and may provide beneficial effects against inflammation and oxidative stress to improve biocompatibility.

The unique fixation device is designed to ensure that Bi-Flow ECMO remains correctly positioned for the duration of the procedure.

BI-FLOW ECMO is validated up to 29 days.

ORDER CODE	303-019
SIZE	19 Fr
OUTER DIAMETER (MM)	6.33
CONNECTION SIZE	3/8"
EFFECTIVE LENGTH (CM)	26,7
AVAILABLE COATED	yes
QTY PER BOX	5
NOTE	Kit Includes: Bi-Flow ECMO 19 Fr, Dilators (12, 19 Fr), Needle 18 Ga, Scalpel, Syringe 5 ml, Guidewire (100 cm, 0.035"), Cannula Introducer, Suture ring, Fixation Device
VALIDATION	29 days