


Advanced Circulatory Support, Mechanical Circulatory Support, ProtekDuo, 16/29Fr, 5140-4629



The ProtekDuo® dual lumen cannula provides the flexibility to do what’s right for each unique patient. It is inserted percutaneously at the right internal jugular vein and it is uniquely positioned to drain deoxygenated blood from the right atrium and return it via the pulmonary artery.

ORDER CODE	5140-4629
SIZE	16/29Fr
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	28 proximal / 46 distal
AVAILABLE COATED	No
QTY PER BOX	1