


Arterial Cannulae, Curved Tip, with Suture Flange, Wire-reinforced Tubing, 14 Fr/ 4.5 mm, A212-45B



Livanova Curved Tip Aortic Arch Cannulae allow for a high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	A212-45B
SIZE	14 Fr/ 4.5 mm
CONNECTOR SIZE	1/4"
EFFECTIVE LENGTH (CM)	23
FEATURE	Connector With Luer Vent
AVAILABLE COATED	Depends on Country Registration
QTY PER BOX	10