


## Arterial Cannulae, Straight Tip, with Suture Collar, Wire-reinforced Tubing, 12 Fr/3.8 mm, A232-38C



Livanova Straight Tip Arterial Cannulae allow for a high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	A232-38C
SIZE	12 Fr/3.8 mm
CONNECTOR SIZE	1/4"
EFFECTIVE LENGTH (CM)	23
FEATURE	Connector Without Vent
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10