


Arterial Cannulae, 135° Curved Tip, with Suture Collar, Wire-reinforced Tubing, 12 Fr/3.8 mm, A272-40N



Livanova 135° Curved Tip Arterial Cannulae allow for a high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|-----------------------|-----------------------------------|
| ORDER CODE | A272-40N |
| SIZE | 12 Fr/3.8 mm |
| CONNECTION SITE | 3/16"-1/4" |
| EFFECTIVE LENGTH (CM) | 22 |
| FEATURE | Without Connector |
| AVAILABLE COATED | Depending on country registration |
| QTY PER BOX | 10 |