



**Arterial Cannulae, Optiflow Arterial, Curved Tip, Non-reinforced Tubing, 21 Fr/ 7.0 mm, A281-70B**



**PERFORMANCE**



Flow (lpm)	Pressure Drop (mmHg) - F21	Pressure Drop (mmHg) - F24
0	0	0
1	~5	~3
2	~10	~6
3	~18	~10
4	~28	~15
5	~40	~22
6	~55	~30
7	~75	~40
8	~105	~55
9	~140	~75

Livanova Flexible Tip Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	A281-70B
SIZE	21 Fr/ 7.0 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	24
FEATURE	Connector with luer vent down
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10