


Arterial Cannulae, Optiflow Arterial, Straight Tip, Wire-reinforced Tubing, 24 Fr/ 8.0 mm, A292-80C



Livanova Optiflow Arterial Cannulae are designed with a three-dimensional basket tip for better flow dispersion and improved hemodynamics. This unique feature provides a native-like blood flow and avoids the so-called jet stream, thus reducing the risk of plaque dislodgement (Assmann et al. Artificial Organs 2014; 39(3): 203-211). Available in 21 and 24 Fr with and without wire-reinforced tubings and optional luer lock for venting. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|-----------------------|-----------------------------------|
| ORDER CODE | A292-80C |
| SIZE | 24 Fr/ 8.0 mm |
| CONNECTOR SIZE | 3/8" |
| EFFECTIVE LENGTH (CM) | 25 |
| FEATURE | Connector Without Vent |
| AVAILABLE COATED | Depending on country registration |
| QTY PER BOX | 10 |