

Cardioplegia Cannulae, Antegrade Aortic Root, with and without Vent Line, 16 Ga (5 Fr), AR-11016



Livanova Antegrade Cardioplegia Aortic Root Cannulae are intended for rapid and secure cardioplegia perfusion of the coronary arteries. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	AR-11016
SIZE	16 Ga (5 Fr)
INSERTION NEEDLE	18 Ga - Pink
EFFECTIVE LENGTH (CM)	15
FEATURE	Without Vent Line
AVAILABLE COATED	no
QTY PER BOX	10