

Cardioplegia Cannulae, Adapters, CA-31010



Livanova Cardioplegia Adapters are designed for use in conjunction with Livanova Antegrade Cardioplegia Aortic Root and Vessel Cannulae for delivery of cardioplegia solution and for heart venting.

ORDER CODE	CA-31010
PRODUCT DESCRIPTION	Multiple Perfusion Adapter, Triplicate with Vessel Cannulae, Cat. No. VC-12000
LENGTH	16.5 in / 41,9 cm
QTY PER BOX	10