

## Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Basket Tip, Stainless Steel Shaft, Female Luer Lock Connector, 12 Fr, CP-12012



Livanova Antegrade Cardioplegia Coronary Ostia Perfusion Cannulae are intended for cardioplegia delivery into the coronary ostia. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	CP-12012
SIZE	12 Fr
CONNECTION SITE	female luer lock
FEATURE	90° Tip
AVAILABLE COATED	no
QTY PER BOX	10