


Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 6.0 mm, CP-21006



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

ORDER CODE	CP-21006
SIZE	6.0 mm
DISTAL CANNULAE TIP SIZE	3.0 mm
FEATURE	Straight Tip
AVAILABLE COATED	no
QTY PER BOX	10