


## Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 6.0 mm, CP-22006



Livanova Antegrade Cardioplegia Coronary Ostia Perfusion Cannulae are intended for cardioplegia delivery into the coronary ostia. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	CP-22006
SIZE	6.0 mm
DISTAL CANNULAE TIP SIZE	3.0 mm
FEATURE	Right AngleTip
AVAILABLE COATED	no
QTY PER BOX	10