

## Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring, Non-reinforced Tubing., 24 Fr/ 8.0 mm, NA-3418



Livanova Flexible Tip Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	NA-3418
SIZE	24 Fr/ 8.0 mm
CONNECTOR SIZE/ACCEPTANCE	3/8"
EFFECTIVE LENGTH (CM)	25
FEATURE	Porous Vent Plug
AVAILABLE COATED	no
QTY PER BOX	10