


Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring, Non-reinforced Tubing., 18 Fr/ 6.0 mm, NA-3426



Livanova Flexible Tip Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|---------------------------|------------------------|
| ORDER CODE | NA-3426 |
| SIZE | 18 Fr/ 6.0 mm |
| CONNECTOR SIZE/ACCEPTANCE | 3/8" |
| EFFECTIVE LENGTH (CM) | 25 |
| FEATURE | Connector without vent |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |