

## Arterial Cannulae, Flexible, Straight Tip with Suture Ring, Non-reinforced Tubing, 24 Fr/ 8.0 mm, NA-4538



Livanova Flexible Tip Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with and without wire reinforced tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	NA-4538
SIZE	24 Fr/ 8.0 mm
CONNECTOR SIZE/ACCEPTANCE	3/8"
EFFECTIVE LENGTH (CM)	25
FEATURE	Connector with luer vent
AVAILABLE COATED	no
QTY PER BOX	10