

## Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, 3D Tip, Stainless Steel Shaft, Female Luer Lock Connector, 15 Fr, P614-50



Livanova Antegrade Cardioplegia Coronary Ostia Perfusion Cannulae are intended for cardioplegia delivery into the coronary ostia. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	P614-50
SIZE	15 Fr
CONNECTION SITE	female luer lock
FEATURE	45° Tip
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10