

## Cardioplegia Cannulae, Antegrade Aortic Root, with Vent Line, with Overpressure Safety Valve, 3.7 mm, R900-06



Livanova Antegrade Cardioplegia Aortic Root Cannulae are intended for rapid and secure cardioplegia perfusion of the coronary arteries. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	R900-06
SIZE	3.7 mm
INSERTION NEEDLE	14 Ga
EFFECTIVE LENGTH (CM)	13
FEATURE	with vent line
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10