


Arterial Cannulae, PureFlex™, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 18 Fr/6.0 mm, RA18-1BP



Livanova PureFlex™ Curved Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with wire reinforced tubing.

ORDER CODE	RA18-1BP
SIZE	18 Fr/6.0 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	21.6
FEATURE	Porous vent plug
AVAILABLE COATED	no
QTY PER BOX	10