

Arterial Cannulae, PureFlex™, Straight Tip, With Suture Collar, Wire-reinforced Tubing, 18 Fr/6.0 mm, RA18-1SC



Livanova PureFlex™ Straight Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with wire reinforced tubing

ORDER CODE	RA18-1SC
SIZE	18 Fr/6.0 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	21.6
FEATURE	Connector without luer vent
AVAILABLE COATED	no
QTY PER BOX	10