

Arterial Cannulae, PureFlex™, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 20 Fr/6.5 mm, RA20-1BP



Livanova PureFlex™ Curved Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with wire reinforced tubing.

ORDER CODE	RA20-1BP
SIZE	20 Fr/6.5 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	21.6
FEATURE	Porous vent plug
AVAILABLE COATED	no
QTY PER BOX	10