

**Arterial Cannulae, PureFlex™, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 22 Fr/ 7.3 mm, RA22-1BC**



Livanova PureFlex™ Curved Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with wire reinforced tubing.

ORDER CODE	RA22-1BC
SIZE	22 Fr/ 7.3 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	21.6
FEATURE	Connector without luer vent
AVAILABLE COATED	no
QTY PER BOX	10