


**Arterial Cannulae, PureFlex™, Curved Tip, With Suture Collar, Wire-reinforced Tubing, 24 Fr/8.0 mm, RA24-1BB**



Livanova PureFlex™ Curved Arterial Cannulae allow for high flow rate and low resistance, as well as easy handling during insertion. They are available with wire reinforced tubing.

ORDER CODE	RA24-1BB
SIZE	24 Fr/8.0 mm
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	21.6
FEATURE	Connector with luer lock
AVAILABLE COATED	no
QTY PER BOX	10