


Cardioplegia Cannulae, Retrograde, Self-Inflating PVC Balloon, 14 Fr, RCS-13314



Livanova Cardioplegia Retrograde Cannulae are intended for cardioplegia delivery into the coronary sinus. They are available with either a self- or manually-inflating balloon.

ORDER CODE	RCS-13314
SIZE	14 Fr
CONNECTOR	Luer lock
EFFECTIVE LENGTH (CM)	30.5
FEATURE	Solid Malleable Stylet with Handle
AVAILABLE COATED	no
QTY PER BOX	5