

Venous Return Cannulae, Triple Stage, 29/29/29 Fr, RTS-11129



Livanova Triple Stage Venous Return Cannulae, with wire-reinforced tubing and with or without connector, feature a third point of drainage for optimal flow performance. The innovative distal tip is softer, yet it adds rigidity to prevent collapse during vacuum assistance procedure.

ORDER CODE	RTS-11129
SIZE	29/29/29 Fr
CONNECTOR SIZE	3/8"
EFFECTIVE LENGTH (CM)	41
FEATURE	With Connector
AVAILABLE COATED	no
QTY PER BOX	10