

Venous Return Cannulae, Triple Stage, 29/37/37 Fr, RTS-13129



Livanova Triple Stage Venous Return Cannulae, with wire-reinforced tubing and with or without connector, feature a third point of drainage for optimal flow performance. The innovative distal tip is softer, yet it adds rigidity to prevent collapse during vacuum assistance procedure.

| | |
|-----------------------|----------------|
| ORDER CODE | RTS-13129 |
| SIZE | 29/37/37 Fr |
| CONNECTOR SIZE | 1/2" |
| EFFECTIVE LENGTH (CM) | 41 |
| FEATURE | With Connector |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |