


## Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 24 Fr, RV-41024



Livanova Single Stage Venous Return Cannulae are provided with straight and right angled tip, with or without reinforced tubing. The wire reinforced tubing mitigates the risk of kinking due to twisting during the procedure, whereas the non-reinforced version is made with kink-resistant, flexible PVC. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	RV-41024
SIZE	24 Fr
TIP LENGTH (CM)	4
CONNECTOR ACCEPTANCE	1/4"-3/8"
AVAILABLE COATED	no
QTY PER BOX	10