


Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 34 Fr, RV-41034



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

ORDER CODE	RV-41034
SIZE	34 Fr
TIP LENGTH (CM)	5,5
CONNECTOR ACCEPTANCE	3/8"
AVAILABLE COATED	no
QTY PER BOX	10