

Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 21.7 Fr, SU-12402



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

| | |
|-----------------------|-----------------|
| ORDER CODE | SU-12402 |
| SIZE | 21.7 Fr |
| CONNECTION SIZE | 1/4" |
| EFFECTIVE LENGTH (CM) | 36,8 |
| FEATURE | Flexible Tubing |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |