

Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 20 Fr, SU-12502



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	SU-12502
SIZE	20 Fr
CONNECTION SIZE	1/4"
EFFECTIVE LENGTH (CM)	36,8
FEATURE	Malleable Tubing
AVAILABLE COATED	no
QTY PER BOX	10