


## Venous Return Cannulae, Single Stage, Straight Bullet Tip, Wire-reinforced Tubing, 34 Fr, V122-34



Livanova Single Stage Venous Return Cannulae are provided with straight and right angled tip, with or without reinforced tubing. The wire reinforced tubing mitigates the risk of kinking due to twisting during the procedure, whereas the non-reinforced version is made with kink-resistant, flexible PVC. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	V122-34
SIZE	34 Fr
CONNECTOR ACCEPTANCE	3/8"
EFFECTIVE LENGTH (CM)	39
FEATURE	Catalog N
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10