


Venous Return Cannulae, Single Stage, Straight Bullet Tip, Wire-reinforced Tubing, 28 Fr, V900-12



Livanova Dual Stage Venous Return Cannulae are designed for venous drainage of the right atrium and vena cava with a single cannula. The wire-reinforced tubing version is a one-piece construction with a flexible drainage basket, and a flexible, yet firm, lighthouse tip. The non-reinforced version is configured with either rigid or flexible drainage tip. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	V900-12
SIZE	28 Fr
CONNECTOR ACCEPTANCE	3/8"
EFFECTIVE LENGTH (CM)	38
FEATURE	With Suture Ring
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10