


## Venous Return Cannulae, Dual Stage, Wire-reinforced Tubing without Obturator and with Connector, 32/40 Fr, V900-49



Livanova Dual Stage Venous Return Cannulae are designed for venous drainage of the right atrium and vena cava with a single cannula. The wire-reinforced tubing version is a one-piece construction with a flexible drainage basket, and a flexible, yet firm, lighthouse tip. The non-reinforced version is configured with either rigid or flexible drainage tip. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	V900-49
SIZE	32/40 Fr
CONNECTOR SIZE	1/2"
EFFECTIVE LENGTH (CM)	41
FEATURE	With Connector
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10