

## Cardioplegia Cannulae, Vessels, VC-11000



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

|                       |                             |
|-----------------------|-----------------------------|
| ORDER CODE            | VC-11000                    |
| EFFECTIVE LENGTH (CM) | 4                           |
| FEATURE               | Without Wings Without Valve |
| AVAILABLE COATED      | no                          |
| QTY PER BOX           | 30                          |