

Cardioplegia Cannulae, Vessels, VC-11010



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

ORDER CODE	VC-11010
EFFECTIVE LENGTH (CM)	6
FEATURE	Without Wings With Valve
AVAILABLE COATED	no
QTY PER BOX	30