


Cardioplegia Cannulae, Vessels, VC-11100



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

ORDER CODE	VC-11100
EFFECTIVE LENGTH (CM)	4
FEATURE	With Wings Without Valve
AVAILABLE COATED	no
QTY PER BOX	30